

## DIA-DUO

Complete Cordless Obturation System



**BUY & GET:**  
1 BOTTLE GP OBTURATORS  
1 DIA-KONDENSOR  
1 DIA-PROSEAL REG. KIT



[4007-2000] Dia-Duo Regular Kit

### Duo-Gun™ Features

- Variable temperature settings: 140°C, 170°C, 200°C
- Temperature reaches 200°C within 15 seconds
- Designed for both left-handed and right-handed individuals
- Operate hand trigger with minimal effort
- Strong and durable tips are available in both 23G (#4007-2120) and 25G (#4007-2121)

### Duo-Pen™ Features

- Tightly compacts and seals all canals incl. lateral canals
- 3 temperature settings: 170°C, 200°C, and 230°C
- Softens, spreads, cuts, and compacts gutta percha
- Quick heating tip reaches 230°C within 1 second
- 5 color-coded tip sizes: #.04-40, .04-45, .04-50, .06-50, .08-55

## DIA ULTRASONIC

Dual-Frequency Ultrasonic Endo Activator

Watch video:



**BUY & GET:**  
1 DIA-ROOT BIO SEALER!  
(VALUE \$75 USD FREE!)



- 1 Dual-frequency: 30 KHz and 45 KHz
- 2 Flexible and durable NiTi tip
- 3 Easy to use and compact
- 4 Fast charging: 3.5 hours for a full charge
- 5 Long lasting: up to 500 cycles
- 6 Powerful and efficient irrigation

[4009-1001] Regular Kit  
Handpiece, charger, tips (2/box), wrench, adapter & cord

## DiaFil™

Light-Cured Nano Hybrid Composite



**BUY 2, GET 1 FREE OR  
BUY 4, GET 4 FREE!**



\*of the same value


- 1 High tensile strength & low shrinkage
- 2 Easy to place/sculpt and does not slump
- 3 Does not stick to instruments
- 4 Adapts well to margins
- 5 Wear & high fracture resistant
- 6 12 shades with remarkable color stability
- 7 Durable restoration w/ long-lasting gloss shine

### SYRINGE REFILL (1 x 4g)

Shades: A1, A2, A3, A3.5, A4, B1, B2, B3, C2, C3, A20, A30

### CAPSULE REFILL (20 x 0.25g)

Shades: A1, A2, A3, A3.5, B2



## Diapex™

Calcium Hydroxide with Iodoform

**BUY 3, GET 1 FREE!**

Formerly Diapex®.  
Different packaging,  
SAME product  
you've trusted and loved:



- 1 Convenient: pre-mixed paste
- 2 Pharmaceutical grade iodoform
- 3 Excellent accessibility to apex
- 4 Antibacterial & bacteriostatic properties
- 5 Highly radiopaque
- 6 Does not harden, Ca(OH)<sub>2</sub> remains active

[A1001-501] Regular Kit 1 x 2g syringe and 20 tips



**DIA-ROOT™ Bio Sealer**  
BioCeramic Root Canal Sealer

**BUY 4, GET 1 FREE!**

- 1 BioCeramic: excellent biocompatibility
- 2 pH 12: highly antibacterial
- 3 Superior bonding to dentin & GP
- 4 Low solubility: zero shrinkage
- 5 Washout resistant: prevent microleakage
- 6 Optimal flowability of 22mm
- 7 Saves time: pre-mixed and pre-filled



[1003-701] Dia-Root Bio Sealer  
2g syringe, 20 tips, 1 cap