

Our Birthday Celebration Continues With the Gift of Savings!

Improve better patient outcomes with longer lasting restorations.



Offers Valid July 1, 2024 - September 30, 2024

All CARBIDES AND SPECIALTY BURS

BUY 4 PACKS GET 1 PACK AT NO CHARGE**

Valid for 5, 10 and 25 Packs and Bulk Packs.

CHOOSE FROM:

Standard Operative Carbides, FG, FGSL, RA, RASL, HP, Great White® Gold, Trimming & Finishing, Bone Cutters, Zekrya, Fissurotomy®, Smart Burs® II

CARBIDE BULK PACK INCLUDING:

FG, RA, HP, GREAT WHITE®, FISSUROTOMY®, TRIMMING AND FINISHING

** Can mix and match shapes to equal or lesser value. Free product to equal or lesser value. Offer excludes Great White® Ultra, EndoGuide® Burs, Endo Safe End and Oral Surgery Burs.



©2024 SS White® All Rights Reserved. SS White®, Fissurotomy® and SmartBurs® II are registered trademarks of SS White Dental. *Market research and test studies available upon request.**No charge goods must be of equal or lesser value. All offers are limited to stock on hand and are not valid on prior purchases. Offers may not be combined with any other sale, promotion, discount, code, coupon and/or offer. Sales promotion choice must be indicated on invoice sent to SS White® on applicable order. SS White® reserves the right to discontinue a promotion at any time. †SS White's 100% guarantee is on all of our products. If you are not completely satisfied with your purchase, SS White® will issue a full refund of the purchase price through your distributor. To claim no charge goods: Mail a copy of your invoice by October 31, 2024 to: SS White Burs, Inc., Attn: Promotions Department, 1145 Towbin Ave., Lakewood, NJ 08701 or email freegoods@sswhitedental.com. Please note that these are end-user doctor promotions and must be redeemed by the qualifying practice. Please allow 4-6 weeks for free delivery. SS White® does its best to ensure that its photos/illustrations are as true to color as possible. However, due to inconsistencies of various monitors, lighting sources, digital photography, and manufacturing lot variations, SS White cannot guarantee that the color that appears on your screen/printed material accurately portrays the true color of the received product. Images are intended as a guide only.

STERILE R Available in Rolls of 25



Better Patient Outcomes
Improved Efficiency
Faster Practice Growth

1145 Towbin Avenue Lakewood, New Jersey 08701
Visit us on the web at www.sswhitedental.com

Unmatched Strength, Unrivalled Cutting Makes SS White® the #1 Choice of Dentists.

Offers Valid July 1, 2024 - September 30, 2024

BUY 4 GET 1 AT NO CHARGE**

**Can mix and match to equal or lesser value • NO CHARGE PRODUCT TO EQUAL OR LESSER VALUE

CHOOSE FROM: BONE CUTTERS, ZEKRYA, T10 AND THE NEW TRIMMING & FINISHING BURS (7283, 48L-010 AND 48L-012)

3C BONE CUTTERS



SS White® Dental Bone Cutters are Strong and Stand Up to the Rigors of Dental Surgery!

- Exceed ISO specifications for neck strength by over 90% reducing the occurrence of breakage*
- Designed for precise trimming and preparation of hard tissue, bone structure, bone lid and apicoectomy procedures
- Sectioning impacted third molars prior to extraction, separating roots, removal of fractured root stump and external sinus lift procedures
- Ultra-smooth vibration free cutting
- The crosscut shape and external flute depth ensures efficient cutting and expels waste material

SS White® Dental
ANNIVERSARY
CELEBRATION **180** YEARS
Better Patient Outcomes | Improved Efficiency | Faster Practice Growth



ZB ZEKRYA BURS



151-Z and 151L-Z (Zekrya) Bur provides exceptionally smooth cutting and greater efficiency for surgical procedures.

Zekrya Bur features a 6 helical bladed cutting edge design and solid carbide construction which delivers exceptional concentricity and reduced vibrations. SS White® Zekrya Burs provide superior cutting performance for oral surgery procedures including sectioning wisdom teeth, separating roots and removing broken root stump. Zekrya Burs can also be used endodontically to reduce cusps on molars to create easy endodontic access.



T10 METAL & CROWN CUTTING BURS



T-10 Metal Cutting Burs for crown removal.

With an improved design of the TriHawk Talon 10 Bur, the SS White® T-10 burs offer greater cutting efficiency and increased durability. The T-10 Bur can cut vertically as well as horizontally, making it ideal for use in procedures where access is limited. Use the T-10 for metal cutting, cavity prep, crown removal as well as creating endodontic access.



TF TRIMMING & FINISHING BURS



12-Blade Finishing Burs

SS White® 12 Blade trimming/finishing bur with premium quality solid tungsten provide faster and smoother cutting. Twelve bladed finishing burs provide a smooth finish on composite, amalgam, enamel and dentin. The flame shape of the 48L and 283 make these burs ideal bur for use in orthodontic debonding procedures. Carbide finishing burs create a smoother surface finish on most dental materials when compared to diamond instruments.



©2024 SS White® All Rights Reserved. SS White® is a registered trademarks of SS White Dental. *Market research and test studies available upon request.**No charge goods must be of equal or lesser value. All offers are limited to stock on hand and are not valid on prior purchases. Offers may not be combined with any other sale, promotion, discount, code, coupon and/or offer. Sales promotion choice must be indicated on invoice sent to SS White® on applicable order. SS White® reserves the right to discontinue a promotion at any time. †SS White's 100% guarantee is on all of our products. If you are not completely satisfied with your purchase, SS White® will issue a full refund of the purchase price through your distributor. To claim no charge goods: Mail a copy of your invoice by October 31, 2024 to: SS White Burs, Inc., Attn: Promotions Department, 145 Towbin Ave., Lakewood, NJ 08701 or email freegoods@sswhitedental.com. Please note that these are end-user doctor promotions and must be redeemed by the qualifying practice. Please allow 4-6 weeks for free delivery. SS White® does its best to ensure that its photos/illustrations are as true to color as possible. However, due to inconsistencies of various monitors, lighting sources, digital photography, and manufacturing lot variations, SS White cannot guarantee that the color that appears on your screen/printed material accurately portrays the true color of the received product. Images are intended as a guide only.